

One-Year Limited Warranty
(One year from date of purchase.)

Upon examination all electronic components if found to be defective in material and or workmanship Will be repaired or replaced at Pacific Pure's option without charge.

The batteries are warranted for 6 months after purchase. (Misuse or accidental damages are not covered).

Exceptions:

An increase in the operational noise level of the motor usually indicates abusive or excessive strain on the scooter which is not covered under this warranty.

Items which may need replacing under normal wear and tear including tires, tubes, belts, bulbs, upholstery, brakes, motor brushes, fuses, etc. Are not covered under this warranty.

Damage caused by abuse, misuse, accident, negligence, improper operation, excessive loading, improper maintenance or storage, or commercial use are not covered.

Any modifications without prior written consent from Pacific Pure will void all warranties.

Return Policy:

All returns are subject to a 15% restocking, handling fee*and must have a return authorization number (A#).

*All shipping cost and/or taxes are non-refundable.

Your local dealer:

Additional Options:



- **Padded Box for rear storage**
(For Model PP-1602 Only)
- **Helmets**



Pacific Pure

offers affordable solutions
in Personal Mobility Scooters.

These Personal Mobility Devices need
No Gas, No Licence and No Insurance,
and they are reasonably priced.

For viewing or a personal demonstration,
contact: **Bill Clark: 705-653-9050**

14 Trent Drive, PO Box 818
Campbellford, Ontario
K0L 1L0
1-800-905-3536

www.PacificPure.ca



MODEL PP-1600

No Gas!
No Licence!
No Insurance!



Battery-Powered Mobility Device

Adjustable seat

Adjustable armrest

Safety non-flip rear wheels

Variable-speed control

Parking brake

Head light with high & low beams

Turn signals with warning beep

Battery life indicator

Speedometer

Left and right hand horn control

Remote alarm

Security wheel lock

Battery-powered mobility scooter

Different colors available: red, silver, black

Wheel size: 16" x 2.5"

Reverse gear

Top speed: <25km/h

Running distance: <50km

Maximum load: 180 kg

Climbing angle: 10°

Battery: 48V/20Ah

Sealed maintenance-free, lead-acid battery

Motor: brushless 500W

Charger: AC 110v Charging time: ≤6h

Product size: 1550*820*1220mm

Gross weight: 75kgs



Lots of storage:
Storage basket
Paper pocket storage
Under seat storage

No Licence or Plates Required...

**Personal Mobility Devices
(Motorized Wheelchairs and
Medical Scooters)**

**Do not require registration, licence
plates, driver's licence or vehicle
insurance.**

**Persons operating motorized
wheelchairs are treated in the same
way as pedestrians.**

**Wheelchairs can be driven by muscular
power or other types of power, and are
designed for and used by people whose
mobility is limited by one or more
conditions or functional impairments.
In general, municipalities establish
by-laws for where wheelchairs can or
cannot be used. Operators should
check with their local municipality to
ensure by-laws permit their use on
sidewalks.**

**A sidewalk should be the first choice for
someone using a wheelchair or medical
scooter. When there is no wheelchair
accessible curb, the person should
return to the sidewalk at the first
available opportunity.**

**If there is no sidewalk available, people
using wheelchairs or personal mobility
devices should travel, like pedestrians,
along the left shoulder of the roadway
facing oncoming traffic.**

An excerpt from the Ontario Ministry of Transportation
Website <http://www.mto.gov.on.ca/english/dandv/vehicle/emerging/index.shtml#personal>